

Product Datasheet

Recombinant Human HSPB7 Protein, His Tag (orb1552248)

Description	Recombinant Human HSPB7 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	HSPB7
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human HSPB7(Met1-Ile170) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	Q9UBY9
Expiration Date	6 months from date of receipt.