

Product Datasheet

Recombinant Human HGPRT/HPRT1 Protein, His Tag (orb1552265)

Description	Recombinant Human HGPRT/HPRT1 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	HPRT1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human HPRT1(Met1-Ala218) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	P00492
Expiration Date	6 months from date of receipt.