

Product Datasheet

Recombinant Human GFAP Protein, His Tag (orb1552327)

| | |
|------------------------|---------------------------------------------------------------------------------------------------------------------|
| Description | Recombinant Human GFAP Protein, His Tag |
| Species/Host | E. coli |
| Reactivity | Human |
| Endotoxins | Please contact our customer service for more information. |
| Target | GFAP |
| Note | For research use only |
| Purity | ≥90% as determined by SDS-PAGE |
| Source | This protein contains the human GFAP(Leu292-Met432) was fused with the N-terminal His Tag and expressed in E. coli. |
| Uniprot ID | P14136 |
| Expiration Date | 6 months from date of receipt. |