

Product Datasheet

Recombinant Human GDNF Protein, His & Fc Tag (orb1552330)

Description	Recombinant Human GDNF Protein, His & Fc Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	GDNF
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human GDNF(Phe20-Ile211) was fused with the C-terminal His Fc Tag and expressed in Mammalian cells.
Uniprot ID	P39905
Expiration Date	6 months from date of receipt.