

Product Datasheet

Recombinant Human GALE Protein, His Tag (orb1552353)

Description	Recombinant Human GALE Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	GALE
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human GALE(Met1-Ala348) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	Q14376
Expiration Date	6 months from date of receipt.