

Product Datasheet

Recombinant Human GABARAPL1 Protein, His Tag (orb1552358)

Description	Recombinant Human GABARAPL1 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	GABARAPL1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human GABARAPL1(Met1-Lys117) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	Q9H0R8
Expiration Date	6 months from date of receipt.