

Product Datasheet

Recombinant Human GABARAP Protein, His Tag & Fc Tag (orb1552360)

Description	Recombinant Human GABARAP Protein, His Tag & Fc Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	GABARAP
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human GABARAP(Met1-Gln116) was fused with the N-terminal His Tag,C-terminal Fc Tag and expressed in E. coli.
Uniprot ID	O95166
Expiration Date	6 months from date of receipt.