

Product Datasheet

Recombinant Human Fractalkine/CX3CL1 Protein, His Tag (orb1552373)

Description	Recombinant Human Fractalkine/CX3CL1 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	CX3CL1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human CX3CL1(Gln25-Arg339) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	P78423
Expiration Date	6 months from date of receipt.