

Product Datasheet

Recombinant Human FIBP Protein, No Tag (orb1552389)

Description	Recombinant Human FIBP Protein, No Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	FIBP
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human FIBP(Met1-Asp357) was fused without Tag and expressed in E. coli.
Uniprot ID	O43427
Expiration Date	6 months from date of receipt.