

Product Datasheet

Recombinant Mouse Fgf-9 Protein, His Tag (orb1552391)

Description	Recombinant Mouse Fgf-9 Protein, His Tag
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	Please contact our customer service for more information.
Target	Fgf9
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Fgf9(Met1-Ser208) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P54130
Expiration Date	6 months from date of receipt.