

Product Datasheet

Recombinant Human FGF-8E Protein, His Tag (orb1552394)

Description	Recombinant Human FGF-8E Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	FGF8
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human FGF8(Gln23-Arg233) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	P55075
Expiration Date	6 months from date of receipt.