

Product Datasheet

Recombinant Mouse FGF-1/Fgf1 Protein, His Tag (orb1552410)

Description	Recombinant Mouse FGF-1/Fgf1 Protein, His Tag
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	Please contact our customer service for more information.
Target	Fgf1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Fgf1(Phe16-Asp155) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P61148
Expiration Date	6 months from date of receipt.