

Product Datasheet

Recombinant Mouse FcRn/Fcgrt Protein, His Tag (orb1552425)

Description	Recombinant Mouse FcRn/Fcgrt Protein, His Tag
Species/Host	Mammal
Reactivity	Mouse
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	Fcgrt
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Fcgrt(Ser22-Val301&Ile21-Met119) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	Q61559
Expiration Date	6 months from date of receipt.