

Product Datasheet

Recombinant Human FCGR3B Protein, His & Fc Tag (orb1552430)

Description	Recombinant Human FCGR3B Protein, His & Fc Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	FCGR3B
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human FCGR3B(Thr20-Gln208) was fused with the C-terminal His-Fc Tag and expressed in Mammalian cells.
Uniprot ID	O75015
Expiration Date	6 months from date of receipt.