

Product Datasheet

Recombinant Human FCER2/CD23 Protein, His Tag (orb1552432)

Description	Recombinant Human FCER2/CD23 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	FCER2
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human FCER2(Asp48-Ser321) was fused with the N-terminus His Tag and expressed in Mammalian cells.
Uniprot ID	P06734
Expiration Date	6 months from date of receipt.