

Product Datasheet

Recombinant Human FABP1 Protein, His Tag (orb1552454)

Description	Recombinant Human FABP1 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	FABP1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human FABP1(Met1-Ile127) was fused with the N-terminus His Tag and expressed in E. coli.
Uniprot ID	Q6FGL7
Expiration Date	6 months from date of receipt.