

Product Datasheet

Recombinant Human Estrogen Receptor/ESR1 Protein, His Tag (orb1552458)

Description	Recombinant Human Estrogen Receptor/ESR1 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	ESR1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human ESR1(Met1-Gln116) was fused with the N-terminus His Tag and expressed in E. coli.
Uniprot ID	P03372
Expiration Date	6 months from date of receipt.