

Product Datasheet

Recombinant Human ESR2 Protein, His Tag (orb1552460)

Description	Recombinant Human ESR2 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	ESR2
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human ESR2(Met1-Ala323) was fused with the N-terminus His Tag and expressed in E. coli.
Uniprot ID	Q92731
Expiration Date	6 months from date of receipt.