

Product Datasheet

Recombinant Human EIF-4E Protein, His Tag (orb1552501)

Description	Recombinant Human EIF-4E Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	EIF4E
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human EIF4E(Met1-Ile118) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	P06730
Expiration Date	6 months from date of receipt.