

Product Datasheet

Recombinant Human EGR1 Protein, His Tag (orb1552504)

Description	Recombinant Human EGR1 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	EGR1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human EGR1(Gln282-Ser433) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	P18146
Expiration Date	6 months from date of receipt.