

## Product Datasheet

### Recombinant Mouse Egf Protein, His Tag (orb1552509)

|                        |                                                                                                                     |
|------------------------|---------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Recombinant Mouse Egf Protein, His Tag                                                                              |
| <b>Species/Host</b>    | E. coli                                                                                                             |
| <b>Reactivity</b>      | Mouse                                                                                                               |
| <b>Endotoxins</b>      | Please contact our customer service for more information.                                                           |
| <b>Target</b>          | Egf                                                                                                                 |
| <b>Note</b>            | For research use only                                                                                               |
| <b>Purity</b>          | ≥90% as determined by SDS-PAGE                                                                                      |
| <b>Source</b>          | This protein contains the mouse Egf(Asn977-Arg1029) was fused with the C-terminal His Tag and expressed in E. coli. |
| <b>Uniprot ID</b>      | <b>P01132</b>                                                                                                       |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                                      |