

Product Datasheet

Recombinant Human DNA Polymerase/POLD2 Protein, His Tag (orb1552535)

Description	Recombinant Human DNA Polymerase/POLD2 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	POLD2
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human POLD2(Met1-Pro469) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P49005
Expiration Date	6 months from date of receipt.