

Product Datasheet

Recombinant Mouse Dkk-1 Protein, 8*His Tag (orb1552542)

Description	Recombinant Mouse Dkk-1 Protein, 8*His Tag
Species/Host	Mammal
Reactivity	Mouse
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	Dkk1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Dkk1(Thr32-His272) was fused with the N-terminal 8- His Tag and expressed in Mammalian cells.
Uniprot ID	O54908
Expiration Date	6 months from date of receipt.