

Product Datasheet

Recombinant Human DKK-1/DKK1 Protein, 8* His Tag (orb1552543)

Description	Recombinant Human DKK-1/DKK1 Protein, 8* His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	DKK1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human DKK1(Thr32-His266) was fused with the N-terminal 8- His Tag and expressed in Mammalian cells.
Uniprot ID	O94907
Expiration Date	6 months from date of receipt.