

## Product Datasheet

# Recombinant Human Cyclophilin C/PPIC Protein, Trx & His Tag (orb1552571)

<b>Description</b>	Recombinant Human Cyclophilin C/PPIC Protein, Trx & His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	PPIC
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human PPIC(Lys31-Asp182) was fused with the N-terminal Trx, His Tag and expressed in E. coli.
<b>Uniprot ID</b>	<a href="#">P45877</a>
<b>Expiration Date</b>	6 months from date of receipt.