

Product Datasheet

Recombinant *Thermobifida fusca* Cutinase Protein, His Tag (orb1552598)

Description	Recombinant <i>Thermobifida fusca</i> Cutinase Protein, His Tag
Species/Host	<i>E. coli</i>
Reactivity	Bacteria
Endotoxins	Please contact our customer service for more information.
Target	Cutinase
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the <i>Thermomonospora fusca</i> cut-2.KW3(Ala1-Phe261) was fused with the C-terminal His Tag and expressed in <i>E. coli</i> .
Uniprot ID	E5BBQ3
Expiration Date	6 months from date of receipt.