

Product Datasheet

Recombinant Human CSNK1G2 Protein, His Tag (orb1552618)

Description	Recombinant Human CSNK1G2 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	CSNK1G2
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human CSNK1G2(Met18-Lys415) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	P78368
Expiration Date	6 months from date of receipt.