

Product Datasheet

Recombinant Human CPNE1 Protein, His Tag & His Tag (orb1552634)

Description	Recombinant Human CPNE1 Protein, His Tag & His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	CPNE1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human CPNE1(Met1-Ala537) was fused with the C-terminal His Tag and N-terminal His Tag and expressed in E. coli.
Uniprot ID	Q99829
Expiration Date	6 months from date of receipt.