

Product Datasheet

Recombinant Mouse Complement Component C3a/C3 Protein, No Tag (orb1552648)

Description	Recombinant Mouse Complement Component C3a/C3 Protein, No Tag
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	Please contact our customer service for more information.
Target	C3
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse C3(Ser671-Arg748) was fused without Tag and expressed in E. coli.
Uniprot ID	P01027
Expiration Date	6 months from date of receipt.