

Product Datasheet

Recombinant Human CDKN1B Protein, His Tag (orb1552709)

Description	Recombinant Human CDKN1B Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	CDKN1B
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human CDKN1B(Met1-Thr198) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	Q6I9V6
Expiration Date	6 months from date of receipt.