

Product Datasheet

Recombinant Human CDK4 Protein, His & T7 Tag (orb1552713)

Description	Recombinant Human CDK4 Protein, His & T7 Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	CDK4
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human CDK4(Met1-Glu303) was fused with the N-terminal His and T7 Tag and expressed in E. coli.
Uniprot ID	P11802
Expiration Date	6 months from date of receipt.