

## Product Datasheet

# Recombinant Human CDK2 Protein, His Tag (orb1552716)

<b>Description</b>	Recombinant Human CDK2 Protein, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	CDK2
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human CDK2(Met1-Leu298) was fused with the N-terminal His Tag and expressed in E. coli.
<b>Uniprot ID</b>	<a href="#">P24941</a>
<b>Expiration Date</b>	6 months from date of receipt.