

Product Datasheet

Recombinant Human CCL5/RANTES Protein, Fc & His Tag (orb1552730)

Description	Recombinant Human CCL5/RANTES Protein, Fc & His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	CCL5
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human CCL5(Ser24-Ser91) was fused with the C-terminal Fc and His Tag and expressed in Mammalian cells.
Uniprot ID	P13501
Expiration Date	6 months from date of receipt.