

Product Datasheet

Recombinant Mouse Ccl3/MIP-1a Protein, His Tag (orb1552732)

Description	Recombinant Mouse Ccl3/MIP-1a Protein, His Tag
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	Please contact our customer service for more information.
Target	Ccl3
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Ccl3(Ala24-Ala92) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	P10855
Expiration Date	6 months from date of receipt.