

Product Datasheet

Recombinant Human CCL3/MIP-1 alpha Protein, No Tag (orb1552733)

Description	Recombinant Human CCL3/MIP-1 alpha Protein, No Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	CCL3
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human CCL3(Ser24-Ala92) was fused without Tag and expressed in E. coli.
Uniprot ID	P10147
Expiration Date	6 months from date of receipt.