

Product Datasheet

Recombinant Human CCL14a/HCC-1/CCL14 Protein, His Tag (orb1552747)

Description	Recombinant Human CCL14a/HCC-1/CCL14 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	CCL14
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human CCL14(Thr20-Asn93) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	Q16627
Expiration Date	6 months from date of receipt.