

Product Datasheet

Recombinant Human Carbonic Anhydrase VIII/CA8 Protein, His Tag (orb1552775)

Description	Recombinant Human Carbonic Anhydrase VIII/CA8 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	CA8
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human CA8(Ala2-Gln290) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P35219
Expiration Date	6 months from date of receipt.