

Product Datasheet

Recombinant Human Carbonic Anhydrase IV/CA4 Protein, His Tag (orb1552778)

Description	Recombinant Human Carbonic Anhydrase IV/CA4 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	CA4
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human CA4(Ala19-Lys283) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P22748
Expiration Date	6 months from date of receipt.