

Product Datasheet

Recombinant Human CAMK1D Protein, Gst Tag (orb1552783)

Description	Recombinant Human CAMK1D Protein, Gst Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	CAMK1D
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human CAMK1D(Met1-Lys385) was fused with the N-terminal Gst Tag and expressed in E. coli.
Uniprot ID	Q8IU85
Expiration Date	6 months from date of receipt.