

Product Datasheet

Recombinant Human Cadherin-9/CDH9 Protein, His Tag (orb1552794)

Description	Recombinant Human Cadherin-9/CDH9 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	CDH9
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human CDH9(Ile22-Ala615) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	Q9ULB4
Expiration Date	6 months from date of receipt.