

## Product Datasheet

# Recombinant Human Betatrophin/ANGPTL8 Protein, His Tag (orb1552841)

<b>Description</b>	Recombinant Human Betatrophin/ANGPTL8 Protein, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	ANGPTL8
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human ANGPTL8(Ala22-Ala198) was fused with the C-terminal His Tag and expressed in E. coli.
<b>Uniprot ID</b>	<a href="#">Q6UXH0</a>
<b>Expiration Date</b>	6 months from date of receipt.