

Product Datasheet

Recombinant Human BCL2L1/BCL-XI Protein, His Tag (orb1552855)

Description	Recombinant Human BCL2L1/BCL-XI Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	BCL2L1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human BCL2L1(Met1-Arg212) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	Q07817
Expiration Date	6 months from date of receipt.