

Product Datasheet

Recombinant Soybean BBI Protein, His Tag (orb1552860)

Description	Recombinant Soybean BBI Protein, His Tag
Species/Host	E. coli
Reactivity	Plant
Endotoxins	Please contact our customer service for more information.
Target	BBI
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the Soybean BBI(Asp40-Asn110) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P01055
Expiration Date	6 months from date of receipt.