

Product Datasheet

Recombinant Rabbit B7-2/CD86 Protein, His Tag (orb1552872)

Description	Recombinant Rabbit B7-2/CD86 Protein, His Tag
Species/Host	Mammal
Reactivity	Rabbit
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	CD86
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the rabbit CD86(Ala23-Arg246) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	P42071
Expiration Date	6 months from date of receipt.