

Product Datasheet

Recombinant Rat B7-1/Cd80 Protein, His Tag (orb1552875)

Description	Recombinant Rat B7-1/Cd80 Protein, His Tag
Species/Host	Mammal
Reactivity	Rat
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	Cd80
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the rat Cd80(Vla30-Gln248) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	G3V671
Expiration Date	6 months from date of receipt.