

Product Datasheet

Recombinant Human ATG3 Protein, T8 Tag & His Tag (orb1552891)

Description	Recombinant Human ATG3 Protein, T8 Tag & His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	ATG3
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human ATG3(Met1-His287) was fused with the N-terminal T8 Tag,C-terminal His Tag and expressed in E. coli.
Uniprot ID	Q9NT62
Expiration Date	6 months from date of receipt.