

Product Datasheet

Recombinant Human AREG Protein, No Tag (orb1552906)

Description	Recombinant Human AREG Protein, No Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	AREG
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human AREG(Ser101-Lys198) was fused without Tag and expressed in E. coli.
Uniprot ID	P15514
Expiration Date	6 months from date of receipt.