

Product Datasheet

Recombinant Human APT-1/LYPLA1 Protein, His Tag (orb1552907)

Description	Recombinant Human APT-1/LYPLA1 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	LYPLA1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human LYPLA1(Met1-Asp230) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	O75608
Expiration Date	6 months from date of receipt.