

## Product Datasheet

# Recombinant Human Activin RIIA/ACVR2A Protein, His Tag (orb1552972)

<b>Description</b>	Recombinant Human Activin RIIA/ACVR2A Protein, His Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1.0 EU per µg of the protein as determined by the LAL method.
<b>Target</b>	ACVR2A
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human ACVR2A(Ala20-Pro1354) was fused with the C-terminal His Tag and expressed in Mammalian cells.
<b>Uniprot ID</b>	<b>P27038</b>
<b>Expiration Date</b>	6 months from date of receipt.