

Product Datasheet

Recombinant Human Activin RIIA/ACVR2A Protein, Fc&His Tag (orb1552973)

Description	Recombinant Human Activin RIIA/ACVR2A Protein, Fc&His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	ACVR2A
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human ACVR2A(Ala20-Pro134) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	P27037
Expiration Date	6 months from date of receipt.